

PagLab Protocols

Laboratory of Dr. Pagliassotti
Department of Food Science and Human Nutrition
Colorado State University
www.paglab.com

TG Analysis of Liver

1. Homogenize approximately 20 mg of liver in 1 ml PBS/10 mM EDTA pH 7.4.
2. Add 2 ml of isopropanol-hexane-water (IHW: 80:20:2, v/v/v) to glass tubes.
 - a. 78.4 : 19.6 : 2.0 ml for 100 ml total volume.
 - b. Glass tubes should be soaked overnight in dd H₂O and dried.
3. Add 200 ul of homogenate to each glass tube from #2 above.
 - a. Mix by vortexing.
 - b. Incubate covered with aluminum foil for 30 minutes at room temperature.
4. Add 500 ul hexane-ether (1:1) to each tube.
 - a. 5 : 5 ml for 10 ml total volume.
 - b. Mix by vortexing.
 - c. Incubate covered with aluminum foil for 10 minutes at room temperature.
5. Add 1 ml of dd H₂O to separate phases.
 - a. Mix by vortexing.
 - b. Incubate covered with aluminum foil for 20 minutes at room temperature.
6. Remove 900 ul of organic phase and place into glass GC vials.
 - a. Evaporate under N₂ gas.
7. Add 500 ul of Infinity TG reagent to each vial, including vials containing TG standard dilutions (0, 0.03125, 0.0625, 0.125, 0.25, 0.5, 1; 12.5 ul standard per 500ul TG reagent)
 - a. Cap vials.
 - b. Incubate at 37°C for 90 minutes with shaking.
 - c. After incubation, dilute 5-fold.
8. Transfer 200 ul from vial to 96-well plate.
9. Determine A₅₄₀ against standard curve (TG standard solution = 200mg/dl).